

Product Datasheet

Human O-6-Methylguanine DNA Methyltransferase (MGMT) ELISA Kit (orb1205229)

Description	Human O-6-Methylguanine DNA Methyltransferase (MGMT) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P16455
Sensitivity	0.051ng/mL
Entrez	4255
Expiration Date	Please enquire.