

Product Datasheet

Human Methyl CpG Binding Protein 2 (MECP2) ELISA Kit (orb1205237)

Description	Human Methyl CpG Binding Protein 2 (MECP2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P51608
Sensitivity	0.065ng/mL
Entrez	4204
Expiration Date	Please enquire.