

Product Datasheet

Human 5-Hydroxytryptamine Receptor 1A (HTR1A) ELISA Kit (orb1205383)

Description	Human 5-Hydroxytryptamine Receptor 1A (HTR1A) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P08908
Sensitivity	0.059ng/mL
Entrez	3350
Expiration Date	Please enquire.