

Product Datasheet

Human High Mobility Group AT Hook Protein 2 (HMGA2) ELISA Kit (orb1205416)

Description	Human High Mobility Group AT Hook Protein 2 (HMGA2) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P52926
Sensitivity	0.114ng/mL
Entrez	8091
Expiration Date	Please enquire.