

Product Datasheet

Human Hemoglobin Gamma 1 (HBg1) ELISA Kit (orb1205434)

Description	Human Hemoglobin Gamma 1 (HBg1) ELISA Kit
Reactivity	Human
Range	468-30000ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or erythrocyte lysates.
Uniprot ID	P69891
Sensitivity	122.6ng/mL
Entrez	3047
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)