

Product Datasheet

Human Hydroxyacyl Coenzyme A Dehydrogenase Beta (HADHb) ELISA Kit (orb1205440)

Description	Human Hydroxyacyl Coenzyme A Dehydrogenase Beta (HADHb) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P55084
Sensitivity	0.053ng/mL
Entrez	3032
Expiration Date	Please enquire.