

Product Datasheet

Human Glutamate Receptor, Ionotropic, AMPA 2 (GRIA2) ELISA Kit (orb1205458)

Description	Human Glutamate Receptor, Ionotropic, AMPA 2 (GRIA2) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P42262
Sensitivity	0.115ng/mL
Entrez	2891
Expiration Date	Please enquire.