

Product Datasheet

Human G Protein Coupled Receptor 109B (GPR109B) ELISA Kit (orb1205468)

Description	Human G Protein Coupled Receptor 109B (GPR109B) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P49019
Sensitivity	0.061ng/mL
Entrez	8843
Expiration Date	Please enquire.