

Product Datasheet

Human Glucagon Like Peptide 1 Receptor (GLP1R) ELISA Kit (orb1205483)

Description	Human Glucagon Like Peptide 1 Receptor (GLP1R) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P43220
Sensitivity	0.118ng/mL
Entrez	2740
Expiration Date	Please enquire.