

Product Datasheet

Human Gamma-Glutamyltransferase 5 (gGT5) ELISA Kit (orb1205494)

Description	Human Gamma-Glutamyltransferase 5 (gGT5) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P36269
Sensitivity	0.55ng/mL
Entrez	2687
Expiration Date	Please enquire.