

Product Datasheet

Human Dual Specificity Phosphatase 6 (DUSP6) ELISA Kit (orb1205620)

Description	Human Dual Specificity Phosphatase 6 (DUSP6) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q16828
Sensitivity	0.122ng/mL
Entrez	1848
Expiration Date	Please enquire.