

Product Datasheet

Human Dual Specificity Phosphatase 1 (DUSP1) ELISA Kit (orb1205622)

Description	Human Dual Specificity Phosphatase 1 (DUSP1) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P28562
Sensitivity	12.9pg/mL
Entrez	1843
Expiration Date	Please enquire.