

## Product Datasheet

# Human Defensin Beta 103A (DEFb103A) ELISA Kit (orb1205668)

<b>Description</b>	Human Defensin Beta 103A (DEFb103A) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	46.8-3000pg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Uniprot ID</b>	<b>P81534</b>
<b>Sensitivity</b>	12.55pg/mL
<b>Entrez</b>	<b>414325</b>
<b>Expiration Date</b>	Please enquire.