

Product Datasheet

Human Defensin Alpha 6, Paneth Cell Specific (DEFa6) ELISA Kit (orb1205669)

Description	Human Defensin Alpha 6, Paneth Cell Specific (DEFa6) ELISA Kit
Reactivity	Human
Range	2.34-150ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q01524
Sensitivity	1.07ng/mL
Entrez	1671
Expiration Date	Please enquire.