

## Product Datasheet

### Human Decorin (DCN) ELISA Kit (orb1205678)

<b>Description</b>	Human Decorin (DCN) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1.25-80ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<b>P07585</b>
<b>Sensitivity</b>	0.43ng/mL
<b>Entrez</b>	<b>1634</b>
<b>Expiration Date</b>	Please enquire.