

Product Datasheet

Human Diazepam Binding Inhibitor (DBI) ELISA Kit (orb1205683)

Description	Human Diazepam Binding Inhibitor (DBI) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P07108
Sensitivity	0.054ng/mL
Entrez	1622
Expiration Date	Please enquire.