

Product Datasheet

Human Corticotropin Releasing Hormone Receptor 2 (CRHR2) ELISA Kit (orb1205743)

Description	Human Corticotropin Releasing Hormone Receptor 2 (CRHR2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q13324
Sensitivity	0.051ng/mL
Entrez	1395
Expiration Date	Please enquire.