

Product Datasheet

Human Cholinergic Receptor, Nicotinic, Alpha 5 (CHRNA5) ELISA Kit (orb1205818)

Description	Human Cholinergic Receptor, Nicotinic, Alpha 5 (CHRNA5) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P30532
Sensitivity	0.048ng/mL
Entrez	1138
Expiration Date	Please enquire.