

Product Datasheet

Human Cholinergic Receptor, Nicotinic, Alpha 3 (CHRNA3) ELISA Kit (orb1205819)

| | |
|--------------------------|---|
| Description | Human Cholinergic Receptor, Nicotinic, Alpha 3 (CHRNA3) ELISA Kit |
| Reactivity | Human |
| Range | 0.312-20ng/mL |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | tissue homogenates or other biological fluids. |
| Uniprot ID | P32297 |
| Sensitivity | 0.107ng/mL |
| Entrez | 1136 |
| Expiration Date | Please enquire. |