

Product Datasheet

Human Calprotectin (CALPRO) ELISA Kit (orb1205870)

Description	Human Calprotectin (CALPRO) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P06702, P05109
Sensitivity	13.12pg/mL
Entrez	6279
Expiration Date	Please enquire.