

Product Datasheet

Human Baculoviral IAP Repeat Containing Protein 7 (BIRC7) ELISA Kit (orb1205896)

Description	Human Baculoviral IAP Repeat Containing Protein 7 (BIRC7) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q96CA5
Sensitivity	0.073ng/mL
Entrez	79444
Expiration Date	Please enquire.