

Product Datasheet

Human A Disintegrin And Metalloprotease 15 (ADAM15) ELISA Kit (orb1206028)

Description	Human A Disintegrin And Metalloprotease 15 (ADAM15) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q13444
Sensitivity	0.054ng/mL
Entrez	8751
Expiration Date	Please enquire.