

## Product Datasheet

# Human ATP Binding Cassette Transporter C8 (ABCC8) ELISA Kit (orb1206039)

<b>Description</b>	Human ATP Binding Cassette Transporter C8 (ABCC8) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.312-20ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	tissue homogenates, cell lysates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q09428</a>
<b>Sensitivity</b>	0.128ng/mL
<b>Entrez</b>	<a href="#">6833</a>
<b>Expiration Date</b>	Please enquire.