

Product Datasheet

Human ATP Binding Cassette Transporter C2 (ABCC2) ELISA Kit (orb1206044)

Description	Human ATP Binding Cassette Transporter C2 (ABCC2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q92887
Sensitivity	0.067ng/mL
Entrez	1244
Expiration Date	Please enquire.