

## Product Datasheet

# Human ATP Binding Cassette Transporter C11 (ABCC11) ELISA Kit (orb1206047)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | Human ATP Binding Cassette Transporter C11 (ABCC11) ELISA Kit |
| <b>Reactivity</b>        | Human   |
| <b>Range</b>             | 0.312-20ng/mL   |
| <b>Note</b>              | For research use only   |
| <b>Application notes</b> | sandwich  |
| <b>Sample Types</b>      | tissue homogenates, cell lysates or other biological fluids.  |
| <b>Uniprot ID</b>        | <a href="#">Q96J66</a>  |
| <b>Sensitivity</b>       | 0.103ng/mL  |
| <b>Entrez</b>            | <a href="#">85320</a>   |
| <b>Expiration Date</b>   | Please enquire.   |