

Product Datasheet

Human ATP Binding Cassette Transporter A4 (ABCA4) ELISA Kit (orb1206048)

Description	Human ATP Binding Cassette Transporter A4 (ABCA4) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P78363
Sensitivity	0.112ng/mL
Entrez	24
Expiration Date	Please enquire.