

## Product Datasheet

# Mouse Transferrin Receptor (TFR) ELISA Kit (orb1206141)

<b>Description</b>	Mouse Transferrin Receptor (TFR) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.781-50ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q62351</a>
<b>Sensitivity</b>	0.28ng/mL
<b>Entrez</b>	<a href="#">22042</a>
<b>Expiration Date</b>	Please enquire.