

Product Datasheet

Mouse Thioredoxin Binding Protein 2 (TBP2) ELISA Kit (orb1206148)

Description	Mouse Thioredoxin Binding Protein 2 (TBP2) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q6SJ95
Sensitivity	0.114ng/mL
Entrez	227606
Expiration Date	Please enquire.