

Product Datasheet

Mouse Surfactant Associated Protein D (SPD) ELISA Kit (orb1206156)

Description	Mouse Surfactant Associated Protein D (SPD) ELISA Kit
Reactivity	Mouse
Range	0.625-40ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, lung lavage fluid, tissue homogenates or other biological fluids.
Uniprot ID	P50404
Sensitivity	0.285ng/mL
Entrez	20390
Expiration Date	Please enquire.