

## Product Datasheet

# Rat Prostaglandin E Synthase, Microsomal (PTGES) ELISA Kit (orb1206208)

<b>Description</b>	Rat Prostaglandin E Synthase, Microsomal (PTGES) ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	0.156-10ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	tissue homogenates, cell lysates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q9JHF3</a>
<b>Sensitivity</b>	0.068ng/mL
<b>Entrez</b>	<a href="#">59103</a>
<b>Expiration Date</b>	Please enquire.