

Product Datasheet

Mouse Progesterone Receptor (PGR) ELISA Kit (orb1206244)

Description	Mouse Progesterone Receptor (PGR) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	A6H6A5
Sensitivity	0.128ng/mL
Entrez	18667
Expiration Date	Please enquire.