

## Product Datasheet

# Rat Cross Linked N-Telopeptide Of Type I Collagen (NTXI) ELISA Kit (orb1206274)

<b>Description</b>	Rat Cross Linked N-Telopeptide Of Type I Collagen (NTXI) ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	3.12-200ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	competitive inhibition
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Sensitivity</b>	0.88ng/mL
<b>Expiration Date</b>	Please enquire.