

## Product Datasheet

# Rat N-Terminal Pro Atrial Natriuretic Peptide (NT-ProANP) ELISA Kit (orb1206277)

<b>Description</b>	Rat N-Terminal Pro Atrial Natriuretic Peptide (NT-ProANP) ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	0.312-20ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Sensitivity</b>	0.148ng/mL
<b>Expiration Date</b>	Please enquire.