

Product Datasheet

Rat Neurotrophin 3 (NT3) ELISA Kit (orb1206281)

Description	Rat Neurotrophin 3 (NT3) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P18280
Sensitivity	0.063ng/mL
Entrez	81737
Expiration Date	Please enquire.