

Product Datasheet

Rat Myelin Oligodendrocyte Glycoprotein (MOG) ELISA Kit (orb1206322)

Description	Rat Myelin Oligodendrocyte Glycoprotein (MOG) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q63345
Sensitivity	0.057ng/mL
Entrez	24558
Expiration Date	Please enquire.