

Product Datasheet

Rat Macrophage Derived Chemokine (MDC) ELISA Kit (orb1206350)

Description	Rat Macrophage Derived Chemokine (MDC) ELISA Kit
Reactivity	Rat
Range	23.4-1500pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q5I0L5
Sensitivity	8.45pg/mL
Entrez	117551
Expiration Date	Please enquire.