

Product Datasheet

Mouse Mannose Associated Serine Protease 1 (MASP1) ELISA Kit (orb1206358)

Description	Mouse Mannose Associated Serine Protease 1 (MASP1) ELISA Kit
Reactivity	Mouse
Range	2.34-150ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	P98064
Sensitivity	0.8ng/mL
Entrez	17174
Expiration Date	Please enquire.