

## Product Datasheet

### Rat Kallikrein 6 (KLK6) ELISA Kit (orb1206391)

<b>Description</b>	Rat Kallikrein 6 (KLK6) ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	0.312-20ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<b>O54854</b>
<b>Sensitivity</b>	0.109ng/mL
<b>Entrez</b>	<b>29245</b>
<b>Expiration Date</b>	Please enquire.