

## Product Datasheet

### Rat Kallikrein 4 (KLK4) ELISA Kit (orb1206392)

<b>Description</b>	Rat Kallikrein 4 (KLK4) ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	0.156-10ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q6IE12</a>
<b>Sensitivity</b>	0.052ng/mL
<b>Entrez</b>	<a href="#">408210</a>
<b>Expiration Date</b>	Please enquire.