

Product Datasheet

Mouse Glycoprotein Ib Beta Polypeptide, Platelet (GP1Bb) ELISA Kit (orb1206455)

Description	Mouse Glycoprotein Ib Beta Polypeptide, Platelet (GP1Bb) ELISA Kit
Reactivity	Mouse
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	plasma.
Uniprot ID	P56400
Sensitivity	14.7pg/mL
Entrez	14724
Expiration Date	Please enquire.