

Product Datasheet

Mouse Gonadotropin Releasing Hormone Receptor (GnRHR) ELISA Kit (orb1206457)

Description	Mouse Gonadotropin Releasing Hormone Receptor (GnRHR) ELISA Kit
Reactivity	Mouse
Range	0.39-25ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q01776
Sensitivity	0.116ng/mL
Entrez	14715
Expiration Date	Please enquire.