

Product Datasheet

Mouse Receptor I For The Fc Region Of Immunoglobulin G (FcγRI) ELISA Kit (orb1206491)

Description	Mouse Receptor I For The Fc Region Of Immunoglobulin G (FcγRI) ELISA Kit
Reactivity	Mouse
Range	0.625-40ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P26151
Sensitivity	0.257ng/mL
Entrez	14129
Expiration Date	Please enquire.