

Product Datasheet

Rat Chemokine C-X-C-Motif Receptor 3 (CXCR3) ELISA Kit (orb1206537)

Description	Rat Chemokine C-X-C-Motif Receptor 3 (CXCR3) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9JII9
Sensitivity	0.059ng/mL
Entrez	84475
Expiration Date	Please enquire.