

Product Datasheet

Rat Cyclin Dependent Kinase 5 (CDK5) ELISA Kit (orb1206580)

Description	Rat Cyclin Dependent Kinase 5 (CDK5) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q03114
Sensitivity	0.059ng/mL
Entrez	140908
Expiration Date	Please enquire.