

## Product Datasheet

# Mouse Amyloid Precursor Protein (APP) ELISA Kit (orb1206658)

<b>Description</b>	Mouse Amyloid Precursor Protein (APP) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.468-30ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">P12023</a>
<b>Sensitivity</b>	0.17ng/mL
<b>Entrez</b>	<a href="#">11820</a>
<b>Expiration Date</b>	Please enquire.