

Product Datasheet

Rat Alcohol Dehydrogenase 1 (ADH1) ELISA Kit (orb1206697)

Description	Rat Alcohol Dehydrogenase 1 (ADH1) ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P06757
Sensitivity	0.59ng/mL
Entrez	24172
Expiration Date	Please enquire.