

## Product Datasheet

# Mouse Alpha-1-Antitrypsin (a1AT) ELISA Kit (orb1206715)

<b>Description</b>	Mouse Alpha-1-Antitrypsin (a1AT) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	31.2-2000ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<b>P22599</b>
<b>Sensitivity</b>	11.5ng/mL
<b>Entrez</b>	<b>20701</b>
<b>Expiration Date</b>	Please enquire.