

## Product Datasheet

# Porcine Heat Shock 70kDa Protein 1A (HSPA1A) ELISA Kit (orb1206722)

<b>Description</b>	Porcine Heat Shock 70kDa Protein 1A (HSPA1A) ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	3.12-200ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q6S4N2</a>
<b>Sensitivity</b>	1.14ng/mL
<b>Entrez</b>	<a href="#">396648</a>
<b>Expiration Date</b>	Please enquire.