

Product Datasheet

Bovine Beta-Lactoglobulin (bLg) ELISA Kit (orb1206748)

Description	Bovine Beta-Lactoglobulin (bLg) ELISA Kit
Reactivity	Bovine
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P02754
Sensitivity	0.073ng/mL
Entrez	280838
Expiration Date	Please enquire.