

Product Datasheet

Human Interleukin 1 Family, Member 9 (IL1F9) ELISA Kit (orb1206751)

Description	Human Interleukin 1 Family, Member 9 (IL1F9) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9NZH8
Sensitivity	5.9pg/mL
Entrez	56300
Expiration Date	Please enquire.